

蛍光試薬	ロングパス	バンドパス
Acridine Orange + DNA	19002SXJ	39002SXJ
Acridine Orange + DNA and RNA	19002SXJ	
Alexa Fluor 350™	19000SXJ	39000SXJ
Alexa Fluor 488™	1900SXJ	39002SXJ
Alexa Fluor 546™	19004SXJ	39004SXJ
Alexa Fluor 555™	19004SXJ	39004SXJ
Alexa Fluor 568™	19006SXJ	39010SXJ
Alexa Fluor 594™	19006SXJ	39010SXJ
Atto 488	19002SXJ	39002SXJ
Atto 550	19004SXJ	39004SXJ
Auramine-rhodamine	19002SXJ	
Azami Green	19002SXJ	39002SXJ
Biosearch Blue	19000SXJ	39000SXJ
BODIPY FL/pH7.2	19002SXJ	39002SXJ
CAL Fluor® Red 590	19004SXJ	39004SXJ
CAL Fluor® Red 610	19006SXJ	39010SXJ
Calcein	19002SXJ	39002SXJ
Calcium Green™-1	19002SXJ	39002SXJ
Calcofluor White	19000SXJ	
Coumarin	19000SXJ	39000SXJ
Cy2™	19002SXJ	39002SXJ
Cy3.5™	19006SXJ	39010SXJ
Cy3™	19004SXJ	39004SXJ
DAPI	19000SXJ	39000SXJ
Dil	19004SXJ	39004SXJ
DiO	19002SXJ	39002SXJ
DsRed	19004SXJ	39004SXJ
DyLight 350	19000SXJ	39000SXJ
DyLight 488	19002SXJ	39002SXJ
DyLight 549	19004SXJ	39004SXJ
DyLight 594	19006SXJ	39010SXJ
EGFP	19002SXJ	39002SXJ
Emerald GFP	19002SXJ	39002SXJ
FAM	19002SXJ	39002SXJ
FITC	19002SXJ	39002SXJ
Fluo-4	19002SXJ	39002SXJ
FM™ 4-64	19006SXJ	
FusionRed	19006SXJ	39010SXJ
GFP	19002SXJ	39002SXJ
Hoechst 33258	19000SXJ	39000SXJ
JC-1 non-ratiometric	19002SXJ	
Killer Red	19006SXJ	39010SXJ
LysoTracker Blue/MeOH	19000SXJ	39000SXJ
LysoTracker Green/pH 5.2	19002SXJ	39002SXJ
mCherry	19006SXJ	39010SXJ
MitoTracker Green FM/MeOH	19002SXJ	39002SXJ
MitoTracker Orange/MeOH	19004SXJ	39004SXJ
MitoTracker Red/MeOH	19006SXJ	39010SXJ
mKate2	19006SXJ	39010SXJ
mPlum	19006SXJ	
mRFP1	19006SXJ	39010SXJ
mWasabi	19002SXJ	39002SXJ
Oregon Green™ 488	19002SXJ	39002SXJ
R-phycoerythrin	19004SXJ	
ReAsH-CCPGCC	19006SXJ	39010SXJ
Resorufin	19004SXJ	39004SXJ
Rhod-2	19004SXJ	39004SXJ
Rhodamine Red™-X	19006SXJ	39010SXJ
ROX	19006SXJ	39010SXJ
Sulforhodamine 101	19006SXJ	39010SXJ
SYBR® Green I	19002SXJ	39002SXJ
SYTO 9/DNA	19002SXJ	39002SXJ
TagRFP	19004SXJ	39004SXJ
TAMRA	19004SXJ	39004SXJ
tdTomato	19004SXJ	39004SXJ
Tetramethylrhodamine isothiocyanate	19004SXJ	39004SXJ
Texas Red®	19006SXJ	39010SXJ
Texas Red®-X	19006SXJ	39010SXJ
TRITC	19004SXJ	39004SXJ
TurboFP650	19006SXJ	
X-rhod-1/Ca2+	19006SXJ	39010SXJ
ZsGreen1	19002SXJ	39002SXJ