

蛍光試薬	ロングパス		バンドパス	
	シングルチャンネルモデル名	マルチチャンネルファイル タキュープ名	シングルチャンネルモデル名	マルチチャンネルモデル 名
Acridine Orange + DNA	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Acridine Orange + DNA and RNA	NX-4100FLS-19002	FLC-19002B		
Alexa Fluor 350™	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
Alexa Fluor 488™	NX-4100FLS-1900	FLC-1900B	NX-4100FLS-39002	FLC-39002B
Alexa Fluor 546™	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Alexa Fluor 555™	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Alexa Fluor 568™	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Alexa Fluor 594™	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Atto 488	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Atto 550	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Auramine-rhodamine	NX-4100FLS-19002	FLC-19002B		
Azami Green	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Biosearch Blue	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
BODIPY FL/pH7.2	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
CAL Fluor® Red 590	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
CAL Fluor® Red 610	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Calcein	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Calcium Green™-1	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Calcofluor White	NX-4100FLS-19000	FLC-19000B		
Coumarin	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
Cy2™	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Cy3.5™	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Cy3™	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
DAPI	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
Dil	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
DiO	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
DsRed	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
DyLight 350	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
DyLight 488	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
DyLight 549	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
DyLight 594	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
EGFP	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Emerald GFP	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
FAM	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
FITC	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Fluo-4	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
FM™ 4-64	NX-4100FLS-19006	FLC-19006B		
FusionRed	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
GFP	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Hoechst 33258	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
JC-1 non-ratiometric	NX-4100FLS-19002	FLC-19002B		
Killer Red	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
LysoTracker Blue/MeOH	NX-4100FLS-19000	FLC-19000B	NX-4100FLS-39000	FLC-39000B
LysoTracker Green/pH 5.2	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
mCherry	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
MitoTracker Green FM/MeOH	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
MitoTracker Orange/MeOH	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
MitoTracker Red/MeOH	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
mKate2	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
mPlum	NX-4100FLS-19006	FLC-19006B		
mRFP1	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
mWasabi	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
Oregon Green™ 488	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
R-phycoerythrin	NX-4100FLS-19004	FLC-19004B		
ReAsH-CCPGCC	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Resorufin	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Rhod-2	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Rhodamine Red™-X	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
ROX	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Sulforhodamine 101	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
SYBR® Green I	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
SYTO 9/DNA	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B
TagRFP	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
TAMRA	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
tdTomato	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Tetramethylrhodamine isothiocyanate	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
Texas Red®	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
Texas Red®-X	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
TRITC	NX-4100FLS-19004	FLC-19004B	NX-4100FLS-39004	FLC-39004B
TurboFP650	NX-4100FLS-19006	FLC-19006B		
X-rhod-1/Ca2+	NX-4100FLS-19006	FLC-19006B	NX-4100FLS-39010	FLC-39010B
ZsGreen1	NX-4100FLS-19002	FLC-19002B	NX-4100FLS-39002	FLC-39002B