

	ロングパス	バンドパス
蛍光試薬	フィルタキューブ名	フィルタキューブ名
Acridine Orange + DNA	FLC-19002A	FLC-39002A
Acridine Orange + DNA and RNA	FLC-19002A	
Alexa Fluor 350™	FLC-19000A	FLC-39000A
Alexa Fluor 488™	FLC-19002A	FLC-39002A
Alexa Fluor 546™	FLC-19004A	FLC-39004A
Alexa Fluor 555™	FLC-19004A	FLC-39004A
Alexa Fluor 568™	FLC-19006A	FLC-39010A
Alexa Fluor 594™	FLC-19006A	FLC-39010A
Atto 488	FLC-19002A	FLC-39002A
Atto 550	FLC-19004A	FLC-39004A
Auramine-rhodamine	FLC-19002A	
Azami Green	FLC-19002A	FLC-39002A
Biosearch Blue	FLC-19000A	FLC-39000A
BODIPY FL/pH7.2	FLC-19002A	FLC-39002A
CAL Fluor® Red 590	FLC-19004A	FLC-39004A
CAL Fluor® Red 610	FLC-19006A	FLC-39010A
Calcein	FLC-19002A	FLC-39002A
Calcium Green™-1	FLC-19002A	FLC-39002A
Calcofluor White	FLC-19000A	
Coumarin	FLC-19000A	FLC-39000A
Cy2™	FLC-19002A	FLC-39002A
Cy3.5™	FLC-19006A	FLC-39010A
Cy3™	FLC-19004A	FLC-39004A
DAPI	FLC-19000A	FLC-39000A
Dil	FLC-19004A	FLC-39004A
DiO	FLC-19002A	FLC-39002A
DsRed	FLC-19004A	FLC-39004A
DyLight 350	FLC-19000A	FLC-39000A
DyLight 488	FLC-19002A	FLC-39002A
DyLight 549	FLC-19004A	FLC-39004A
DyLight 594	FLC-19006A	FLC-39010A
EGFP	FLC-19002A	FLC-39002A
Emerald GFP	FLC-19002A	FLC-39002A
FAM	FLC-19002A	FLC-39002A
FITC	FLC-19002A	FLC-39002A
Fluo-4	FLC-19002A	FLC-39002A
FM™ 4-64	FLC-19006A	
FusionRed	FLC-19006A	FLC-39010A
GFP	FLC-19002A	FLC-39002A
Hoechst 33258	FLC-19000A	FLC-39000A
JC-1 non-ratiometric	FLC-19002A	
Killer Red	FLC-19006A	FLC-39010A
LysoTracker Blue/MeOH	FLC-19000A	FLC-39000A
LysoTracker Green/pH 5.2	FLC-19002A	FLC-39002A
mCherry	FLC-19006A	FLC-39010A
MitoTracker Green FM/MeOH	FLC-19002A	FLC-39002A
MitoTracker Orange/MeOH	FLC-19004A	FLC-39004A
MitoTracker Red/MeOH	FLC-19006A	FLC-39010A
mKate2	FLC-19006A	FLC-39010A
mPlum	FLC-19006A	
mRFP1	FLC-19006A	FLC-39010A
mWasabi	FLC-19002A	FLC-39002A
Oregon Green™ 488	FLC-19002A	FLC-39002A
R-phycoerythrin	FLC-19004A	
ReAsH-CCPGCC	FLC-19006A	FLC-39010A
Resorufin	FLC-19004A	FLC-39004A
Rhod-2	FLC-19004A	FLC-39004A
Rhodamine Red™-X	FLC-19006A	FLC-39010A
ROX	FLC-19006A	FLC-39010A
Sulforhodamine 101	FLC-19006A	FLC-39010A
SYBR® Green I	FLC-19002A	FLC-39002A
SYTO 9/DNA	FLC-19002A	FLC-39002A
TagRFP	FLC-19004A	FLC-39004A
TAMRA	FLC-19004A	FLC-39004A
tdTomato	FLC-19004A	FLC-39004A
Tetramethylrhodamine isothiocyanate	FLC-19004A	FLC-39004A
Texas Red®	FLC-19006A	FLC-39010A
Texas Red®-X	FLC-19006A	FLC-39010A
TRITC	FLC-19004A	FLC-39004A
TurboFP650	FLC-19006A	
X-rhod-1/Ca2+	FLC-19006A	FLC-39010A
ZsGreen1	FLC-19002A	FLC-39002A