

蛍光試薬	ロングパス	バンドパス
	シングルチャンネルモデル名	シングルチャンネルモデル名
Acridine Orange + DNA	NX-4100FLS-19002	NX-4100FLS-39002
Acridine Orange + DNA and RNA	NX-4100FLS-19002	
Alexa Fluor 350™	NX-4100FLS-19000	NX-4100FLS-39000
Alexa Fluor 488™	NX-4100FLS-19002	NX-4100FLS-39002
Alexa Fluor 546™	NX-4100FLS-19004	NX-4100FLS-39004
Alexa Fluor 555™	NX-4100FLS-19004	NX-4100FLS-39004
Alexa Fluor 568™	NX-4100FLS-19006	NX-4100FLS-39010
Alexa Fluor 594™	NX-4100FLS-19006	NX-4100FLS-39010
Atto 488	NX-4100FLS-19002	NX-4100FLS-39002
Atto 550	NX-4100FLS-19004	NX-4100FLS-39004
Auramine-rhodamine	NX-4100FLS-19002	
Azami Green	NX-4100FLS-19002	NX-4100FLS-39002
Biosearch Blue	NX-4100FLS-19000	NX-4100FLS-39000
BODIPY FL/pH7.2	NX-4100FLS-19002	NX-4100FLS-39002
CAL Fluor® Red 590	NX-4100FLS-19004	NX-4100FLS-39004
CAL Fluor® Red 610	NX-4100FLS-19006	NX-4100FLS-39010
Calcein	NX-4100FLS-19002	NX-4100FLS-39002
Calcium Green™-1	NX-4100FLS-19002	NX-4100FLS-39002
Calcofluor White	NX-4100FLS-19000	
Coumarin	NX-4100FLS-19000	NX-4100FLS-39000
Cy2™	NX-4100FLS-19002	NX-4100FLS-39002
Cy3.5™	NX-4100FLS-19006	NX-4100FLS-39010
Cy3™	NX-4100FLS-19004	NX-4100FLS-39004
DAPI	NX-4100FLS-19000	NX-4100FLS-39000
Dil	NX-4100FLS-19004	NX-4100FLS-39004
DiO	NX-4100FLS-19002	NX-4100FLS-39002
DsRed	NX-4100FLS-19004	NX-4100FLS-39004
DyLight 350	NX-4100FLS-19000	NX-4100FLS-39000
DyLight 488	NX-4100FLS-19002	NX-4100FLS-39002
DyLight 549	NX-4100FLS-19004	NX-4100FLS-39004
DyLight 594	NX-4100FLS-19006	NX-4100FLS-39010
EGFP	NX-4100FLS-19002	NX-4100FLS-39002
Emerald GFP	NX-4100FLS-19002	NX-4100FLS-39002
FAM	NX-4100FLS-19002	NX-4100FLS-39002
FITC	NX-4100FLS-19002	NX-4100FLS-39002
Fluo-4	NX-4100FLS-19002	NX-4100FLS-39002
FM™ 4-64	NX-4100FLS-19006	
FusionRed	NX-4100FLS-19006	NX-4100FLS-39010
GFP	NX-4100FLS-19002	NX-4100FLS-39002
Hoechst 33258	NX-4100FLS-19000	NX-4100FLS-39000
JC-1 non-ratiometric	NX-4100FLS-19002	
Killer Red	NX-4100FLS-19006	NX-4100FLS-39010
LysoTracker Blue/MeOH	NX-4100FLS-19000	NX-4100FLS-39000
LysoTracker Green/pH 5.2	NX-4100FLS-19002	NX-4100FLS-39002
mCherry	NX-4100FLS-19006	NX-4100FLS-39010
MitoTracker Green FM/MeOH	NX-4100FLS-19002	NX-4100FLS-39002
MitoTracker Orange/MeOH	NX-4100FLS-19004	NX-4100FLS-39004
MitoTracker Red/MeOH	NX-4100FLS-19006	NX-4100FLS-39010
mKate2	NX-4100FLS-19006	NX-4100FLS-39010
mPlum	NX-4100FLS-19006	
mRFP1	NX-4100FLS-19006	NX-4100FLS-39010
mWasabi	NX-4100FLS-19002	NX-4100FLS-39002
Oregon Green™ 488	NX-4100FLS-19002	NX-4100FLS-39002
R-phycoerythrin	NX-4100FLS-19004	
ReAsH-CCPGCC	NX-4100FLS-19006	NX-4100FLS-39010
Resorufin	NX-4100FLS-19004	NX-4100FLS-39004
Rhod-2	NX-4100FLS-19004	NX-4100FLS-39004
Rhodamine Red™-X	NX-4100FLS-19006	NX-4100FLS-39010
ROX	NX-4100FLS-19006	NX-4100FLS-39010
Sulforhodamine 101	NX-4100FLS-19006	NX-4100FLS-39010
SYBR® Green I	NX-4100FLS-19002	NX-4100FLS-39002
SYTO 9/DNA	NX-4100FLS-19002	NX-4100FLS-39002
TagRFP	NX-4100FLS-19004	NX-4100FLS-39004
TAMRA	NX-4100FLS-19004	NX-4100FLS-39004
tdTomato	NX-4100FLS-19004	NX-4100FLS-39004
Tetramethylrhodamine isothiocyanate	NX-4100FLS-19004	NX-4100FLS-39004
Texas Red®	NX-4100FLS-19006	NX-4100FLS-39010
Texas Red®-X	NX-4100FLS-19006	NX-4100FLS-39010
TRITC	NX-4100FLS-19004	NX-4100FLS-39004
TurboFP650	NX-4100FLS-19006	
X-rhod-1/Ca2+	NX-4100FLS-19006	NX-4100FLS-39010
ZsGreen1	NX-4100FLS-19002	NX-4100FLS-39002