

蛍光試薬	ロングパス	バンドパス
	落射蛍光ユニット名	落射蛍光ユニット名
Acridine Orange + DNA	FIS-SXJ-19002	FIS-SXJ-39002
Acridine Orange + DNA and RNA	FIS-SXJ-19002	
Alexa Fluor 350™	FIS-SXJ-19000	FIS-SXJ-39000
Alexa Fluor 488™	FIS-SXJ-19002	FIS-SXJ-39002
Alexa Fluor 546™	FIS-SXJ-19004	FIS-SXJ-39004
Alexa Fluor 555™	FIS-SXJ-19004	FIS-SXJ-39004
Alexa Fluor 568™	FIS-SXJ-19006	FIS-SXJ-39010
Alexa Fluor 594™	FIS-SXJ-19006	FIS-SXJ-39010
Atto 488	FIS-SXJ-19002	FIS-SXJ-39002
Atto 550	FIS-SXJ-19004	FIS-SXJ-39004
Auramine-rhodamine	FIS-SXJ-19002	
Azami Green	FIS-SXJ-19002	FIS-SXJ-39002
Biosearch Blue	FIS-SXJ-19000	FIS-SXJ-39000
BODIPY FL/pH7.2	FIS-SXJ-19002	FIS-SXJ-39002
CAL Fluor® Red 590	FIS-SXJ-19004	FIS-SXJ-39004
CAL Fluor® Red 610	FIS-SXJ-19006	FIS-SXJ-39010
Calcein	FIS-SXJ-19002	FIS-SXJ-39002
Calcium Green™-1	FIS-SXJ-19002	FIS-SXJ-39002
Calcofluor White	FIS-SXJ-19000	
Coumarin	FIS-SXJ-19000	FIS-SXJ-39000
Cy2™	FIS-SXJ-19002	FIS-SXJ-39002
Cy3.5™	FIS-SXJ-19006	FIS-SXJ-39010
Cy3™	FIS-SXJ-19004	FIS-SXJ-39004
DAPI	FIS-SXJ-19000	FIS-SXJ-39000
Dil	FIS-SXJ-19004	FIS-SXJ-39004
DiO	FIS-SXJ-19002	FIS-SXJ-39002
DsRed	FIS-SXJ-19004	FIS-SXJ-39004
DyLight 350	FIS-SXJ-19000	FIS-SXJ-39000
DyLight 488	FIS-SXJ-19002	FIS-SXJ-39002
DyLight 549	FIS-SXJ-19004	FIS-SXJ-39004
DyLight 594	FIS-SXJ-19006	FIS-SXJ-39010
EGFP	FIS-SXJ-19002	FIS-SXJ-39002
Emerald GFP	FIS-SXJ-19002	FIS-SXJ-39002
FAM	FIS-SXJ-19002	FIS-SXJ-39002
FITC	FIS-SXJ-19002	FIS-SXJ-39002
Fluo-4	FIS-SXJ-19002	FIS-SXJ-39002
FM™ 4-64	FIS-SXJ-19006	
FusionRed	FIS-SXJ-19006	FIS-SXJ-39010
GFP	FIS-SXJ-19002	FIS-SXJ-39002
Hoechst 33258	FIS-SXJ-19000	FIS-SXJ-39000
JC-1 non-ratiometric	FIS-SXJ-19002	
Killer Red	FIS-SXJ-19006	FIS-SXJ-39010
LysoTracker Blue/MeOH	FIS-SXJ-19000	FIS-SXJ-39000
LysoTracker Green/pH 5.2	FIS-SXJ-19002	FIS-SXJ-39002
mCherry	FIS-SXJ-19006	FIS-SXJ-39010
MitoTracker Green FM/MeOH	FIS-SXJ-19002	FIS-SXJ-39002
MitoTracker Orange/MeOH	FIS-SXJ-19004	FIS-SXJ-39004
MitoTracker Red/MeOH	FIS-SXJ-19006	FIS-SXJ-39010
mKate2	FIS-SXJ-19006	FIS-SXJ-39010
mPlum	FIS-SXJ-19006	
mRFP1	FIS-SXJ-19006	FIS-SXJ-39010
mWasabi	FIS-SXJ-19002	FIS-SXJ-39002
Oregon Green™ 488	FIS-SXJ-19002	FIS-SXJ-39002
R-phycoerythrin	FIS-SXJ-19004	
ReAsH-CCPGCC	FIS-SXJ-19006	FIS-SXJ-39010
Resorufin	FIS-SXJ-19004	FIS-SXJ-39004
Rhod-2	FIS-SXJ-19004	FIS-SXJ-39004
Rhodamine Red™-X	FIS-SXJ-19006	FIS-SXJ-39010
ROX	FIS-SXJ-19006	FIS-SXJ-39010
Sulforhodamine 101	FIS-SXJ-19006	FIS-SXJ-39010
SYBR® Green I	FIS-SXJ-19002	FIS-SXJ-39002
SYTO 9/DNA	FIS-SXJ-19002	FIS-SXJ-39002
TagRFP	FIS-SXJ-19004	FIS-SXJ-39004
TAMRA	FIS-SXJ-19004	FIS-SXJ-39004
tdTomato	FIS-SXJ-19004	FIS-SXJ-39004
Tetramethylrhodamine isothiocyanate	FIS-SXJ-19004	FIS-SXJ-39004
Texas Red®	FIS-SXJ-19006	FIS-SXJ-39010
Texas Red®-X	FIS-SXJ-19006	FIS-SXJ-39010
TRITC	FIS-SXJ-19004	FIS-SXJ-39004
TurboFP650	FIS-SXJ-19006	
X-rhod-1/Ca2+	FIS-SXJ-19006	FIS-SXJ-39010
ZsGreen1	FIS-SXJ-19002	FIS-SXJ-39002